FACS Protocol: Surface and Intracellular Staining of Human Whole Blood

- A. Surface staining
 - 4 Add 100 μl of whole blood in 5ml BD falcon (FACS tube).
 - Add 10 μl of each antibody (commercial) in every sample. Add appropriate isotypic controls.
 - **4** Mix gently and incubate in a dark place for *15 min* at room temperature.
 - **4** Add 2ml **BD lysing solution 1x** in every sample.
 - ↓ Vortex for 10 sec.
 - **4** Incubate in a dark place for *15 min* at room temperature.
 - **4** Centrifuge at 1000 rpm for *5 min*.
 - Discard supernatants.
 - 4 Add 1ml **PBS 1x** without Ca or Mg.
- B. Intracellular staining
 - For each FACS sample use 0.5 2 x 10⁶ cells/tube.
 Wash cells by adding 3 ml 5% PBS/FCS per sample.
 Centrifuge at 200 x g for 5 min.
 - 4 Aspiration (100-150 μl) Discard supernatants.
 - Fix by adding 100 μl Reagent A (monoclonal-mouse)
 Mix gently by hand and incubate for *15 min* at room temperature.
 - Wash cells by adding 3 ml 5% PBS/FCS.
 Centrifuge at 200 x g for 5 min.
 - **4** Aspiration Discard supernatant leaving some covering the cells.
 - 4 Permeabilization Add 100 μl Reagent B and mix gently.
 - Add 1.8 μl of Antibody for TF.
 Spin
 Incubate for 3 h at 4 °C.
 - Wash cells by adding 3 ml PBS/FCS and mix gently. Centrifuge at 200 x g for 5 min.
 - ♣ Aspiration Discard supernatants.

Add 10 μl secondary Antibody (anti-mouse) and incubate for 30 min at 4°C. Vortex.

Wash cells by adding 3ml 5% PBS/FCS.

Centrifuge at 200 x g for 5 min.